eQo 1-Step ToughMix®

No compromises. Unparalleled performance meets sustainability

Shipping

FEATURES AND BENEFITS:

- Unrivaled Sensitivity Detection down to 2 copies of RNA
- Sustainable & Earth Friendly lyophilized eQo bead format eliminates dry ice, CO, consumptions & reduces shipping costs
- Novel Reverse Transcriptase convenient, ambient reaction setup. Overcomes rigid RNA secondary structure
- High Stability Up to 1 year at room temperature, up to 1 year rehydrated at –20°C
- Tough Tested Broad and significantly enhanced tolerance to PCR inhibitors

DESCRIPTION:

eQo 1-step ToughMix is a lyophilized reagent system for reverse transcription quantitative PCR (RT-qPCR) of RNA templates using hybridization probe detection chemistries such as TaqMan® 5'-hydrolysis probes. It is supplied with a proprietary rehydration buffer that when combined with the lyophilized "eQo beads" produces a stabilized 4X concentrated master mix. The kit includes thermolabile UDG for amplicon carryover elimination, an enhanced warm-start reverse transcriptase (RT) and all other required components for 1-step RT-qPCR except RNA template and primer/probe(s). The reaction chemistry has been optimized for inhibitor tolerance and delivers exceptional performance in either single-plex or highly demanding multiplex 1-step RT-qPCR formats.

Enhanced Thermostability

The included qScript Ultra reverse transcriptase supports rapid and processive first-strand synthesis at temperatures up to 65°C (optimal 55°C to 60°C), which disrupts interfering RNA secondary structure and improves primer specificity. This novel RT is further enhanced by an aptamer "warm-start" component

that effectively blocks RT activity during reaction setup enabling highly sensitive and reproducible low copy quantification and extended room-temperature stability of fully assembled reactions.






Figure 1 Amplification of a virioid RNA with a high degree of stable secondary structure. Dahlia latent virioid RNA template was prepared, and upon folding, the ssRNA molecule was shown to adopt a stable rod-like structure. Compared with kits from other suppliers, reactions with eQo 1-step carried out at a reverse transcription temperature of 60°C show consistently improved detection.






Overcome Challenging Samples

eQo 1-Step ToughMix demonstrates a tolerance to a broad range of common PCR inhibitors, reducing the requirement for perfectly "clean" samples. Additionally, amplification of difficult gene targets, include those that are GC-rich, is robust and consistent.





Inhibitor resistance across wide range of sample types. Innate inhibitor tolerance of qScript Ultra reverse transcriptase together with the formulation of the eQo 1-step mix allow for superior tolerance to many common reaction inhibitors. Reverse transcription was carried out at manufacturer's recommendations (including 60°C for eQo 1-Step); thermocycling was performed at the same cycling protocols for all. Cq values were compared with reactions without inhibitors to show relative inhibitor tolerance in the 1-step reaction.

Figure 3 Gene expression quantification using primer limiting of GC-rich targets. eQo 1-step shows consistent mRNA quantification over a range of difficult gene targets. 1-step reaction mixtures were assembled in duplicate using 50 ng total human RNA and the indicated primer/probe assay sets. Reverse transcription was carried out at manufacturer's recommendations (including 60°C for eQo 1-Step); thermocycling was performed at the same cycling protocols for all The average Cq was plotted for each gene expression assay.

Extreme Sensitivity

Detection of extremely low-copy RNA virus is possible with eQo 1-Step ToughMix. Amplification of just 2 copies of SARS-CoV-2 RNA template is robust and clearly evident.




Figure 4 Improved limits of detection for RNA pathogens: In a fourplex RT-qPCR experiment, eQo 1-step ToughMix can detect RNA pathogens in a purified sample matrix with as few as two copies. A representative amplification plot was generated by diluting SARS-CoV-2 RNA from 200,000 to 2 copies. The inset plot displays the average Cq against the log copy number of the SARS-CoV-2 target. Strong amplification was observed even at the highest dilution.

ORDER INFO

Product Name	Quantabio Catalog Number	Size
eQo 1-Step ToughMix - 100 R	95301-100	100 x 20 μL rxns
eQo 1-Step ToughMix - 1000 R	95301-02K	1000 x 20 μL rxns

