

**NEST**<sup>®</sup>

Tissue Culture &  
Lab Plasticware

**CELL CULTURE**  
Dishes, Plates & Flasks  
cell scrapers  
**CENTRIFUGE TUBES**  
serological & pasteur pipettes  
**CRYOGENIC VIALS**  
Glass Bottomed  
**CELL CULTURE DISHES**  
**Reagent Reservoirs**  
**DEEP WELL PLATES**



Report No.01893



ISO9001:2008



Registration No. 10035061





ISO 13485:2003

## Cell Culture Dishes



Cat #	Description	Pack size
706001	35mm Diameter cell culture dish, sterile, inner pack of 20 x 25	500
705001	60mm Diameter cell culture dish, sterile, inner pack of 20 x 25	500
704001	100mm Diameter cell culture dish, sterile, inner pack of 20 x 15	300

## Cell Culture Plates



Cat #	Description	Pack size
703001	6-well Cell culture plate, flat, clear, individually wrapped, sterile	50
712001	12-well Cell culture plate, flat, clear, individually wrapped, sterile	50
702001	24-well Cell culture plate, flat, clear, individually wrapped, sterile	50
748001	48-well Cell culture plate, flat, clear, individually wrapped, sterile	50
701001	96-well Cell culture plate, flat, clear, individually wrapped, sterile	100

## Cell Culture Flasks



### Plug Seal Cap

Cat #	Description	Pack size
707001	25cm <sup>2</sup> Cell culture flask, 50ml volume, sterile, inner pack of 10 x 20	200
708001	75cm <sup>2</sup> Cell culture flask, 250ml volume, sterile, inner pack of 5 x 20	100
709001	175cm <sup>2</sup> Cell culture flask, 750ml volume, sterile, inner pack of 5 x 8	40

### Filter Vent Cap

Cat #	Description	Pack size
707003	25cm <sup>2</sup> Cell culture flask, 50ml volume, sterile, inner pack of 10 x 20	200
708003	75cm <sup>2</sup> Cell culture flask, 250ml volume, sterile, inner pack of 5 x 20	100
709003	175cm <sup>2</sup> Cell culture flask, 750ml volume, sterile, inner pack of 5 x 8	40

## Cell Scrapers



Cat #	Description	Pack size
710001	Cell scraper, handle 220mm, free rotating blade 13mm, sterile	100
710011	Cell scraper, handle 280mm, free rotating blade 20mm, sterile	100

## Centrifuge Tubes



Cat #	Description	Pack size
601002	15ml Centrifuge Tubes bulk, sterile, bags of 50 x 10	500
601001	15ml Centrifuge tubes in racks of 50 x 10, sterile	500
602002	50ml Centrifuge Tubes bulk, sterile, bags of 25 x 20	500
602001	50ml Centrifuge tubes in racks of 25 x 20, sterile	500
602072	50ml Self standing centrifuge tubes, sterile, resealable bags of 25 x 20	500

## Cryogenic Vials



Cat #	Description	Pack size
618001	0.5ml External thread, round bottom, self-standing, sterile	500
618101	0.5ml Internal thread, round bottom, self-standing, sterile	500
606001	1.2ml External thread, round bottom, self-standing, sterile	500
606101	1.2ml Internal thread, round bottom, self-standing, sterile	500
607001	2.0ml External thread, round bottom, self-standing, sterile	500
607101	2.0ml Internal thread, round bottom, self-standing, sterile	500
608001	4.0ml External thread, round bottom, self-standing, sterile	250
608101	4.0ml Internal thread, round bottom, Self-Standing, sterile	250
609001	5.0ml External thread, round bottom, self-standing, sterile	250
609101	5.0ml Internal thread, round bottom, self-standing, sterile	250

\* Variety of coloured cap inserts available

## Serological Pipettes



Cat #	Description	Pack size
324001	1ml Serological pipette, 280mm, yellow, individually wrapped, sterile	500
325001	2ml Serological pipette, 280mm, green, individually wrapped, sterile	400
326001	5ml Serological pipette, 340mm, blue, individually wrapped, sterile - NEW pack size	200
327001	10ml Serological pipette, 339mm, orange, individually wrapped, sterile	200
328001	25ml Serological pipette, 308.5mm, red, individually wrapped, sterile	200
329001	50ml Serological pipette, 343.5mm, purple, individually wrapped, sterile	100

## Pasteur Pipettes



Cat #	Description	Pack size
318012	1ml Pasteur pipette, 160mm, individually wrapped	500
318212	3ml Pasteur pipette, 162mm, individually wrapped	500

## Deep Well Plates



Cat #	Description	Pack size
501601	96 well, round well plate, 0.5ml well volume, sterile	50
502601	96 well, round well plate, 1ml well volume, sterile	50
503501	96 well, square well plate, 2ml well volume, sterile	50

## Other Products



Cat #	Description	Pack Size
360002	60ml, Reagent reservoir with lid, V bottom	25
801002	35mm Glass bottomed cell culture dishes, 15mm dimension, TC treated, sterile	20
801001	35mm Glass bottomed cell culture dishes, 20mm dimension, TC treated, sterile	20