

MORTEXER™

MORE Power, MORE Flexibility... MORE Vortexer!

- + The only vortexer available with Multi-Head™
- + Powerful Motor for INSTANT Vortexing
- + Q-Drive™ Dynamic Balance System
- + Touch or continuous operation



Unique Multi-Head allows for standard vortexing PLUS holds up to 8 microtubes

The *Mortexer™* is MORE than your standard vortex mixer. The unique design of the Multi-Head™ (included) features a traditional cup head for general purpose mixing, plus holds up to eight microcentrifuge tubes.

The drive system is balanced with Q-Drive™ technology, providing a smooth vortexing action, while minimizing noise and excessive vibration. This results in quieter, more efficient and longer lasting operation. Increased speed, up to 3400 rpm produces an instant vortex of even the largest sample sizes, including full 50ml tubes.

Other features include: “no-walk” suction feet, cold room/incubator compatibility and spill-proof electronic controls. Two year warranty.

Technical Data:

| | |
|------------------------|---|
| Speed Range: | 200 to 3400 rpm |
| Operating Modes: | Touch or continuous |
| Operating Temp. Range: | +4°C to +65°C |
| Dimensions: | 13 x 16 x 17 cm 5 x 6.3 x 6.7 in. |
| Weight: | 3.8 kg / 8.4 lb |
| Electrical: | 115V AC, 60 Hz, 1.5A 230V AC, 50 Hz, 0.75A 230V AC, 60 Hz, 0.75A |
| Warranty: | 2 years |



Ordering Information:

BV1005* *Mortexer™* Vortex Mixer (with Multi-Head)

Optional Accessories:

| | |
|---------------------|---|
| BV1000-FLAT | Optional 3 inch flat head |
| BV1000-COMBO | Optional head for tubes and microplate. 66 tubes (38 x 1.5ml and 28 x 0.5ml) or 1 microplate |
| BV1000-H15 | Optional head for 12 x 1.5/2.0ml, horizontal |
| BV1000-H150 | Optional head for 4 x 15ml, horizontal |
| BV1000-H500 | Optional head for 2 x 50ml, horizontal |

* To order a product in 230V / 50Hz, please add -E to item number.
To order a product in 230V / 60Hz, please add -E2 to item number.